




Explorers Double End

Used for caries and calculus detection
Exploration & Restorations
Handle No. 2 (6mm) or No. 4 (8mm)



Explorers




Probes



Expors



04-1145




Explorers and Probes


A Double end instrument with an Explorer on one end and Probe on the other for charting
Handle No. 2 (6mm) or No. 4 (8mm)



Spreaders

Used to compact filling material during lateral condensation
Handle No.2 (6mm)





04-1181
U of W
Palmer



04-1182
Brewer




04-1183
Palmer




04-1184
Ivory



04-1185
R.D Punch
Alnsworth



04-1189
R.D Punch
Ivory



R.D Frames

Pedo 3 1/2"
04-1186

Child 4"
04-1187

Adult 5"
04-1188



Steiglitz STR
04-1190



Steiglitz 45"
04-1191




Steiglitz 90"
04-1192




Pluggers Anterior


Used to compact filling material during vertical condensation
the instrument are marked at 5mm intervals
Handle No. 2 (6mm)



Pluggers

Used to Compact filling material during
vertical condensation
Handle No. 2 (6mm)







Pluggers Posterior
Used to compact filling material during vertical condensation
The instruments are marked at 5mm intervals
Handle No. 2 (6mm)



Blue Titanium Coated Composite Instruments


Light Weight Silicon Handles with Blue Titanium Coated
Non Stick Tips Use For Composite Placement & Contouring





Blue Titanium Coated Composite Instruments


Light Weight Silicon Handles with Blue Titanium Coated
Non Stick Tips Use For Composite Placement & Contouring







Gold Titanium Coated Composite Instruments

Soft Grip Handles with Gold Titanium Non Stick
Tips Use For Composite Placement & Contouring
Handle No.4 (8mm)





Gold Titanium Coated Burnishers Soft Grip
Soft Grip Handles with Gold Titanium Non Stick
Tips Use For Composite Placement & Contouring
Handle No.4 (8mm)




Gold Titanium Coated Composite Instruments
Soft Grip Handles with Gold Titanium Non Stick
Tips Use For Composite Placement & Contouring
Handle No.4 (8mm)




Hatchets


Use For cavity Preparation, Internal Line Angles
Retentive areas and removing hard caries
Handle No.4 (8mm)



Margin Trimmers and Angle Formers

Use To Produce a Proper Bevel on Enamel Margins
Handle No.4 (8mm)






Excavators English Patterns
For removal of carious dentin
carving amalgam and direct wax patterns
Handle No.4 (8mm)


Excavators Oval Spoons and Blades
For removal of carious dentin,
surretaw tceip dne megleme gniarec rof also



Composite Plastic Filling Instruments

Highly Polished Instruments Tips used
For Composite Placement and Counting





Excavators Spoons

For removal of carious dentin
also for carving amalgam and direct wax patterns
Handle No.4 (8mm)

Composite Plastic Filling Instruments


Highly Polished Instruments Tips used
For Composite Placement and counting
Handle No.4 (8mm)



Carvers

Used to Carve anatomical features and trim excess materials
Handle No.4 (8mm)






Carvers
Used to Carve anatomical
features and trim excess materials
Handle No.4 (8mm)




Carving Knives

Used for trimming Excess Filling material
Flash and overhangs
Handle No.4 (8mm)



Burnishers


Designed to condense, smooth
carve and polish amalgam
Handle No.4 (8mm)






Pluggers Condensers

Used to condense filling materials into cavity preparations
Handle No. 4 (8mm)



Placement Instruments and Gingival Retractors


Placement Instruments Used to Place Base Re-liner with in cavity preparations
Gingival Retractors Protects tissue during cavity preparations and placement of composite restorations
Handle No. 4 (8mm)





Pluggers Condensers

Used to condense filling materials into cavity preparations
Handle No. 4 (8mm)





Straight

04-1960

Ochsenbein

04-1959

Ochsenbein

04-1958

Binangle (15-10-10)

04-1957

Binangle (15-8-8)

04-1956

Binangle (20-9-8)

04-1955

Straight (15-20)

04-1954

Wedelstaedt (15-15-3)


04-1953

Wedelstaedt (11-5-16-3)

04-1952

Wedelstaedt (115-16-3)

04-1951



Amalgam well stainless steel
with non slip base ring

Chisels
Used for refine the cavity preparation
Handle No. 4 (8mm)

Amalgam Carrier
To carry and dispense filling materials
Amalgam, Composite and Bone.



Ortho Instruments

Heavy duty hard wire cutters, face bending pliers and intraoral photographic mirrors



**Side Cutter Tc
Side Cutter
Upto 0.040"**
04-2001



**End Cutter Tc
Side Cutter
Upto 0.040"**
04-2002



**Derer Facebow Plier
Facebow loop forming
in either inner or outer bows
upto 0.063**
04-2003



**Magnum
Universal Facebow
look forming and cutting
plier upto 0.071"**
04-2004



**Inter dental Stripping
Manual Stripper handle only holds all types of abrasive strips**
04-2005



**Heavy duty Cutter
hard wire cutter compound
action stainless steel
upto 0.057"**
04-2006



Lingual
04-2007



Buccal
04-2008



Palatal Adult
04-2009



Palatal Child
04-2010



Loop Forming Schwarz Mini Plier
Optical style. round and concave beaks loop forming plier upto 0.020"
04-2021



Loop Forming Schwarz Midi Plier
Optical style. round and concave beaks loop forming plier upto 0.028"
04-2022



Loop Forming Schwarz Mixi Plier
Optical style. round and concave beaks loop forming plier upto 0.028"
04-2023



Universal Dolphin Style
Used for bending and cutting pliers upto 0.036"
04-2024



Waldsachs
Used for bending and cutting pliers upto 0.036"
04-2025



Clasp Forming Nance Plier
Adams clasp forming plier upto 0.028"
04-2026



Clasp Forming Schwarz Straight Adams for arrow class forming plier upto 0.028"
04-2027



Class Forming Schwarz Angular used for arrow class forming plier upto 0.028"
04-2028



Crimping Plier Straight
used to crimp surgical ball hooks to archwire, slotted tip fits cover base of hook for secure hold while crimping
04-2029



Crimping Plier angled
used to crimp surgical ball hooks to archwire, slotted tip fits cover base of hook for secure hold while crimping
04-2030



Ortho Instruments

Wire bending and loop forming pliers



Young Plier
used for loop
forming
04-2031



Adams Plier mini
used for prices right
angle bends and
appliance
adaptation up to
0.040"
04-2032



Adams Plier midi
used for prices right
angle bends and
appliance
adaptation up to
0.040"
04-2033



Adams Plier mixi
used for prices right
angle bends and
appliance
adaptation up to
0.040"
04-2034



Adams Plier Long
used for prices right
angle bends and
appliance
adaptation up to
0.040"
04-2035



Veingart Fine Tips
Slim, tapered 20° tips
and versatility
04-2036



Veingart Fine Tips TC
Slim, tapered 20° tips
and versatility
04-2037



Veingart Standard Tips
Versatile utility plier tapered
32° tips convenient for
placement and removal of
archwires
04-2038



Veingart Standard Tips TC
Versatile utility plier tapered
32° tips convenient for
placement and removal of
archwires
04-2039



Jarabak Plier
use for light wire bending
with serrated tips and groves
assure accurate bending and
closing of loops
04-2040



**Long Nose Light Wire
Bird Beak Plier**
use for light wire bending
and loops forming plier
long and fine beaks for
making smaller loops
04-2041



**Standard Angle
Bird Plier**
use for light wire bending
and loops forming the tips are
extremely fine for making
small loops
04-2042



**Angle Bird Break Plier
with TC Cutter**
Universal with all wire
bending and loop forming
features of bird beak pliers
plus cutter
04-2043



**Aderer Long Tips
Clasp Adjusting Long
3 Prong Plier for wire
bending & Contouring**
04-2044



Aderer Mini
Clasp Adjusting with
round 3 Prong nose
Plier for wire
bending & Contouring
04-2045



Aderer Midi
Clasp Adjusting with
round 3 Prong nose
Plier for wire
bending & Contouring
04-2046



Aderer Maxi
Clasp Adjusting with
round 3 Prong nose
Plier for wire
bending & Contouring
04-2047



Round Nose Short Beak
used for wire bending
and loop forming
04-2048



Round Nose Long Beak
used for wire bending
and loop forming
04-2049



Ortho Instruments

Arch forming, torquing, contouring and how utility pliers



Lingual Palatal Arch Forming Plier
use to form double-back or triple-back bends in either 0.30" or 0.036" wire for insertion in lingual sheath
04-2050



Tweed/Angle Arch Forming Plier
use to form rectangular arch forming & torquing 1.78mm wide tips
04-2051



Tweed/Angle Arch Forming Plier Slim
use to form rectangular arch forming & torquing 0.64mm wide tips
04-2052



Dela Rose Arch Forming & Contouring Plier
use to form contouring & forming archwires without scoring or torquing
04-2053



Turrets
used for forming archwires without torque, for wires 0.16" - 0.022"
04-2054



How Plier Str.
u0.30mm to 0.50mm
04-2055



How Plier Cvd.
u0.30mm to 0.50mm
04-2056



How Plier Str. TC Tips.
u0.30mm to 0.50mm
04-2057



How Plier Cvd. TC Tips.
u0.30mm to 0.50mm
04-2058



Distal End Cutter with Safety hold.
Minihead with TC insert, long beak and long straight handle cuts 0.30mm to 0.50mm
04-2059



Soft Wire Cutter & Ligature Cutter, 45° angled with TC tips.
Cuts Wire 0.25 mm to 0.38mm
04-2060



Distal End Cutter with Safety hold.
Minihead with TC insert, long beak and long straight handle cuts 0.30mm to 0.50mm
04-2061



Soft Wire Cutter & Ligature Cutter, Standard with TC tips.
Cuts Wire 0.25 mm to 0.38mm
04-2062



Hard Wire Cutter with TC tips.
Cuts Wire 0.30 mm to 0.50mm
04-2063



Hard Wire Cutter 15° angled with TC tips.
Cuts Wire 0.30 mm to 0.50mm
04-2064



Soft Wire Cutter & Ligature Cutter, Micro with TC tips.
Cuts Wire 0.23 mm to 0.30mm
04-2065



Soft Wire Cutter & Ligature Cutter, Mini with TC tips.
Cuts Wire 0.25 mm to 0.38mm
04-2066



Ortho Instruments

Elastic ligature applicators, placing forceps and separators placing pliers
Buckle tube, Distal end bending & Cutting instruments



Separators Placing Pliers
Straight Beak
04-2067



Separators Placing Pliers
Angled Beak
04-2068



Buccal tube convertible cap
removing plier cuts off the
temporary cap on a convertible
first molar buccal tube when
banding of second molar is
required.
04-2069

Spare tube for buccal tube



Distal end cutter standard
cuts wire max 0.30 to 0.50mm
04-2070



Distal end cutter standard
with safety hold TC Insert
cuts wire max 0.30 to 0.50mm
04-2071



Buccal Tube Tweezer
04-2072



**Band Removing
Plier Left.**
04-2073



**Band Removing
Plier Straight.**
04-2074

**Band Removing Plier removes
posterior bands safety and quickly**



**04-2075
Contouring Plier and
adapting bands or crowns**



**04-2076
Contouring Plier and
adapting bands or crowns**



**04-2077
Bracket Removing Plier
Straight.**



**04-2078
Bracket removing Plier 60° angled tip
allows easy access to posterior brackets.**

**Designed curved tips,
provide easy access
and grip of bands
without access retraction
of cheek or lips.**



**04-2079
Bracket Removing Plier removes O.B, Soft
pad rest on incisal surface of tooth, sharp
tip wedge the brackets base.**




**Band Removing
Plier Right.**
**04-2080
Contouring Plier and
adapting bands or crowns**



Ortho Instruments


Measuring Instruments and accessories for brackets



Surgery Knives

Used for incisions that remove soft and interproximal tissue
Handle No. 6 (10mm)






Surgical Chisels
Used for removing and reshaping bones
Handle No. 4 (8mm)



Periosteals

For reflecting and retracting the mucoperiosteum after incisions of the gingival tissue







Universal Curettes
Blade with two cutting edges and rounded toe good for all surfaces
Handle No. 6 (10mm) or Handle No. 4 (8mm)



Sickle Scaler

Designed for the removal of supragingival Calculus
Handle No.6 (10mm) or Handle No. 4 (8mm)





Gracy Curettes Standard
Designed to adapt to a specific area to allow for deep scaling or tooth surface
Handle No.6 (10mm) also available in Handle No. 4 (8mm)


Gracy Curettes With Long Shank Tips
Designed to adapt to a specific area to allow for deep scaling or tooth surface
Handle No.6 (10mm) or Handle No. 4 (8mm)

Terminal Shank is 3mm Longer than Standard Gracy Curette to reach deep in to a pocket




Goldman Fox Scalers

Used for removing heavy supragingival calculus
handle No.4 (8mm) or No.6 (10mm)



Gracy Curettes Short Shank Rigid Tips

Designed to adapt to a specific area or tooth surface
handle No.4 (8mm) or No.6 (10mm)






Small cassette for 5 Instruments
3"x8"x1.25"
76x203x32 mm
04-2441



04-2442
Large cassette for
16 Instruments with an accessory Area
8"x11"x1.25"
203x280x32 mm




04-2445
Large cassette for
20 Instruments
8"x11"x1.25"
203x280x32 mm




04-2443
1 Tissue Forcep 130mm
1 Crile haemostatic forceps 140mm
1 straight surgical scissors 130mm
1 Dissecting forceps 130 mm
1 Feilchenfeld splinter forceps 115mm
25 Michel clip 16*3
1 flacon of suture silk, dterile in alcohol



04-2444
Large cassette for
Elevators & Luxators
8"x11"x2"
203x280x50.8 mm



Automatic Crown Remover
with Tips and 2 Tools
04-2471



Manual Crown Remover
with 3 Tips
04-2477



Contrangle
Junior
04-2479



Lvory 1
04-2478



Universal
04-2480



Contrangle
Senior
04-2481

Instruments Cassettes
Available in Red, Blue, Yellow and Green Silicon
inserts, Also in stainless steel inserts

Crown Removers and Matrix Retainers Bands
Used To remove Crowns



Miscellaneous Utensils

Used to remove crowns



Mallet Hammer with
Stainless Steel Head
7.5 oz
04-2486



Mallet Hammer with
Nylon Caps
8 oz
04-2487



Burs Holder
Mesh
04-2488



Cotton Pellet Dispenser
Stainless Steel
04-2489



Cotton Pellet Dispenser
Plastic
04-2490



Endo Box
04-2491



Kidney Trays
S,M,L,XL
04-2492



Instruments Box
S,M,L
04-2493



Large
04-2494



Medium
04-2495



Small
04-2496



Instruments Box
with handle on lid
S,M,L
04-2497



Upper



04-2501
XL



04-2502
L



04-2503
M



04-2504
S

Lower



XL
04-2505



L
04-2506



M
04-2507



S
04-2508

Upper



XL
04-2509



L
04-2510



M
04-2511



S
04-2512

Lower



XL
04-2513



L
04-2514



M
04-2515



S
04-2516

Impression Trays
Used for taking Impressions of the mouth



Impression Trays
Used for taking Impressions of the mouth

Upper

Solid Trays set of 10



XL
04-2517



L
04-2518



M
04-2519



S
04-2520



XS
04-2521

Lower



XL
04-2522



L
04-2523



M
04-2524



S
04-2525



XS
04-2526

Perforated trays set of 10

Upper



XL
04-2527



L
04-2528



M
04-2529



S
04-2530



XS
04-2531

Lower



XL
04-2532



L
04-2533



M
04-2534



S
04-2535



XS
04-2536



Solid Trays set of 6



U-L
04-2537



U-M
04-2538



U-S
04-2539



L-L
04-2543



L-M
04-2544



L-S
04-2545

Perforated trays set of 6

Upper



U-L
04-2540




U-M
04-2541



U-S
04-2542

Lower



L-L
04-2546



L-M
04-2547



L-S
04-2548



30 Left



31 Right



30 Left



31 Right



32 Lower
04-2549




33 Upper
04-2550



32 Lower
04-2551



33 Upper
04-2552



99 Lateral
04-2553



99 Lateral
04-2554

Perforated set of 5

Solid set of 5


Used for taking Impressions of the mouth

Impression Trays



P.K. Thomas Waxing Carvers

Useful for functional Laboratory waxing





Plaster Nipper
04-2601



Fahanstock Knife Small
04-2602



Fahanstock Knife Large
04-2603



Grtman-10 Plaster Knife,
04-2604



7R Plaster Knife, 2" Blade
04-2605



6R Plaster Knife, 2" Blade
04-2606



15R Rigid, 6" Tapered Blade
04-2607

Plaster Knife & Plaster Nipper
For General Laboratory use

Spatulas



04-2626
24 Flexible
44 mm



04-2627
324 Rigid
51 mm



04-2628
A6 Rigid
25 mm



04-2629
5 Novatech
Long/Fluted



04-2630
6 Novatech
Long/Blade




04-2631
7 Novatech
Short/Blade



04-2632
8 Novatech
Long



04-2633
9 Novatech
Short



04-2636
Spatula Single End
Fig1 Fig2 Fig3 Fig4 Fig5 Fig6



04-2634
7 Wax



04-2635
Spoon / Spatula

Spatulas Wax Knives
Use for functional Laboratory waxing



Mixing Plaster Spatulas

Used for mixing materials and general laboratory uses



1A Stiff,
4 1/2 " Cvd Blade
04-2651



11R Flexible
3 1/2 " Blade
04-2652



10R Extra Stiff
4 1/4 " Blade
04-2653



8R Stiff
4" Blade
04-2654



4R Rigid
4 1/4 " Blade
04-2655



3R Stiff
4" Blade
04-2656



2R Flexible
4" Blade
04-2657

Scalpel Blades Handles and Caliper

Used for cutting, Trimming & Measuring



04-2676



04-2677



04-2678




04-2679



04-2680



04-2681



Boley gauge
04-2682



04-2684


Scalpel Blade No.11

Scalpel Blade No.12

04-2685

Scalpel Blade No.15

04-2686



Lawson Spring Caliper Metal
04-2683



Mixing Cement Spatulas

Used to mix load cement into crowns
or inlay preparations



04-2701

7



04-2702

31



04-2703

Novatech-5



04-2704

Double End



04-2705

324 Rigid



04-2706

24 Flexible




04-2707

22



04-2708

A-6




04-2731

Kramer Nevins Acute




04-2732

Kramer Nevins Obtuse




04-2733

300 Sinus Lift




04-2734

Sinus Lift




04-2735

Sinus Lift



04-2736

1 Sinus Lift 3.5 mm



04-2737

1 Sinus Curette



04-2738

2 Sinus Curette

Sinus Lift Instruments



Osteotomes and Sinus Retractors

Osteotomes



Sinus Retractor
Misch 15mm



Sinus Retractor
Misch 10mm



Sinus Retractor



Osteotome
Javanovic 6mm



Osteotome
Curved 7.5mm



Osteotome
Straight 7.5mm



Osteotome
Curved 4mm



Osteotome
Straight 4mm

Tapered Convex Tip 7-10-13-15-18mm Markings



Osteotome
2.7mm



Osteotome
3.2mm



Osteotome
3.7mm



Osteotome
4.2mm



Osteotome
2.7mm



Osteotome
3.2mm



Osteotome
3.7mm



Osteotome
4.2mm

Tapered Concave Tip 7-10-13-15-18mm Markings

Osteotomes Set



04-2784



04-2785



04-2786



04-2787




04-2788



Bone Scoop 5mm 04-2801




LaBanca 3x1.5mm



04-2802

4 Javanovic 3mm



04-2803

Bone Packer 5/7mm




04-2804

Bone Packer 2/5mm




04-2805

Bone Camper 7mm



04-2806

Bone Rasp



04-2807




04-2808
Bone crusher
Titanium Head




04-2809
Bone Caliper Forceps

colour Probe 20mm




04-2821

Periodontal 3mm/4mm




04-2822

Periodontal 4mm/5mm



04-2823

Periodontal 3mm/5mm




04-2824

Jovanovic 1mm/4mm




04-2825

Jovanovic 2mm/5mm




04-2826

Jovanovic 3mm/6mm



04-2827

Palacci



04-2828

VOLKMANN 21cm



04-2810

UTRINE SCOPE 240 mm




04-2811

Bone Pluggers and Bone Scoop


Periodontal Chisels Implant Chisels and probe




Implant Periosteals



Periotomes



Implant Chisels





Chisel

Single-Bevel
Straight
Scale 9-11-13-15

Single-Bevel
Curved


Double-Bevel
Slightly Curved







Molars General



Upper Molars



Lower Molars



Premolars



Anterior



Order Form

Name		Telephone () - - - -	
Street			
City	Province	Postal Code	- - - - -

For Office Use Only

No	Art No	Description	Qty	Price (CAD\$)	Total

Special Instructions:

Sub Total

GST

Shipping

Total Amount

Payment Method ☐ Credit Card ☐ Cheque or Other

Credit Card Payment Authorization	Total Amount (In Words)		(Please indicate the type of Credit Card)	
			☐ VISA	☐ MasterCard
			☐ AMERICAN EXPRESS	
				
				
Name of the cardholder		Credit Card Number		
Signature of the cardholder		Expiry Date		
		Month Year		
		Date Year Month Day		